



佛山市柯诺医疗器械有限公司
KASO MEDICAL TECHNOLOGY CO.,LTD

地址：佛山市南海区狮山镇罗村沙坑东华工业区
ADD: Donghua Industry Park Luocun Shishan Town, Nanhai District, Foshan City,
Guangdong Province, China

0757-88792672 / 400-0578-810

Tel: +86-757-86361716
Whatsapp: (+86)18666357977
Email: info@kasogroup.com
Facebook: KASO MEDICAL

www.kasomedical.cn

本公司保留不予以事先通知而进行产品改良之权利
Our Company Reserves The Right To Make Any Improvement On Product Design Without Advice Notice


图片仅供参考,以实物为准
All Products Subject To The Practicality, the Pictures Just For Reference



KEEP
ACHIEVING
STAY
OUTSTANDING



Professional Dental Manufacturer



佛山市柯诺医疗器械有限公司
KASO MEDICAL TECHNOLOGY CO.LTD





Kingship Accompanied by **Supreme Optimization**

至尊体验与极致优化于一体

瞰思未来

ENLIGHTENED BY THE FUTURE




KASO Medical Technology Co.,Ltd is a High-tech professional manufacturer specializes in producing precise dental equipment which has the large-scale production facilities,covering a factory area of more than 10000m² and setting manufacturing,sales and after-sales service in one of professional dental comprehensive apparatus manufacturers. KASO also has Ateco RVG sensor imports from UK,and It is the sole agency of Ateco sensor in the world taking charge of marketing and sales in different countries.

KASO has total 19 main departments,including the professional design&development department,new technology training department,quality supervision department,electronic component department,hardware spare part production department,dental unit installation department,X-Ray unit installation department,air compressor installation department,camera LED light department,intraoral camera installation department,inspection department,package department,HR department,legal department,patent department,advertising department,sales&merchandising department,B2B network marketing department and after-sales service department.Every department is in charged by the professional staff who make sure the quality through the strict inspection and by the relative departments which begin from the paper design to the packing.Ourwell-equipped facilities and excellent quality control throughout all stages of production enables us to pass CE,ISO13485,SGS and FSC certificate etc.

In addition,our company regularly participates in World-renowned dental exhibitions such as IDS in Cologne Germany,AEEDC in Dubai UAE,Vietnam International Dental Exhibition,GNYDM in New York USA,Dental-Expo in Moscow Russia, APDC in Manila Philippines,ADEX in Istanbul Turkey, Dental South China exhibition in Guangzhou,SINO-Dental exhibition in Beijing and Dentech China exhibition in Shanghai.We look forward to meeting you at the World Dental Trade Fair to jointly seek

Our company's target is DOING QUALITY SERIOUSLY,MAKING PRODUCTS IRREPLACEABLE. KASO is specialized in manufacturing the high quality and competitive price products through summarizing new opportunities for cooperation,the similar products' strength to meet the market.We take the technique supportive team as the main element in the sells to bring the extreme function of technique to our customers.It has been accepted and supported by lots of customers.KASO will be the leading position in the dental field offering new technology,new ideas and new design products.Therefore,our company's goal is ALWAYS BE IMITATED,NEVER BE SURPASSED.

We have rich export experience and products are sold to 178 countries all over the world,receiving high praise from customers at home and abroad.We look forward to becoming a strategic partner with you.We sincerely hope that we will establish friendly business cooperated relationship with friends come from the global dental field to develop the dental industry.




我们的目标是认真做质量,做到产品无可代替。柯诺通过总结类似产品的实力以满足市场,专门制造高质量和具有竞争力的价格产品。我们以技术支持团队作为销售的主要元素,为客户带来最先进的功能和技术。柯诺已被许多客户接受和支持。KASO 将在牙科领域处于领先地位,提供新技术,新理念和新设计产品。因此,我们的目标是一直被模仿,从未被超越。

我们拥有丰富的出口经验,产品销往全球 178 个国家,受到国内外客户的高度赞誉。我们期待着成为您的战略合作伙伴。我们真诚希望能与来自全球牙科领域的朋友建立友好的业务合作关系,共同发展牙科产业。

佛山市柯诺医疗器械有限公司是一家专业生产精密牙科设备的高新技术专业生产厂家,拥有大规模的生产设备,工厂面积超过 10000 平方米,并将制造、销售和售后服务集于一体的,专业的牙科综合仪器制造商。柯诺还拥有从英国进口的 Ateco RVG 传感器,柯诺是 Ateco 传感器的全球总代理并且负责不同国家市场营销和销售。

柯诺共有 19 个主要部门,包括专业设计研发部门.新技术培训部门.质量监督部门.电子元件部门、硬件配件生产部门、牙科综合治疗台安装部门、X 光机组装部门、空压机安装部门、LED 灯带摄像功能生产部门、口腔内窥镜安装部门、检测部门、包装部门、人力资源部门、法律部门、专利部门、广告部门、销售和跟单部门、B2B 网络营销部门和售后服务部门。每个部门都由专业人员负责,这些专业人员通过严格的检查和相关部门从纸张设计到包装来确保质量。我们在生产的各个阶段都拥有完善的设备和出色的质量控制,使我们能够通过 CE,ISO13485,SGS 和 FSC 等证书。

此外,我司定期参加德国科隆牙科展、迪拜国际口腔医学展览会、美国纽约国际口腔医学大会世界牙科联盟展会、越南国际口腔展、中国广州华南口腔展、北京国际口腔展、上海国际口腔器材展览会、菲律宾亚太牙科展、土耳其国际口腔展览会和俄罗斯牙科博览会等享誉世界的牙科展览会。我们期待在世界牙科贸易展览会上与您会面,共同寻求合作新机遇。



About KASO dental

PROFESSIONAL DENTAL MANUFACTURER

DIGITIZING FUTURE KASO

数字化未来



Global originality 全球独创

- 牙科治疗椅应急供电系统
 - 牙科治疗椅节流变频直流系统
 - 牙科治疗椅纳米净水 UV 消毒系统
 - 工具盘内置消毒盒
- Emergency power supply system for dental chair
 - Restricted-frequency DC system for dental chair
 - Dental chair nano water purification UV disinfection system
 - Tool tray built-in disinfection box



AND STRENGTH WITNESSES MAGNIFICENCE



Business philosophy 经营理念

我们坚信“总有一款适合您！”选择柯诺，
选择您的贴心专业助理

We firmly believe that "there is always one for you!"
Choose KASO, choose your intimate professional
assistant

Enterprise mission 企业使命

立足牙科健康行业，服务全球牙科专家，成为引
领全球的牙科设备综合制造商

Based on the dental health industry, serving the
global dental experts, to become the world's leading
dental equipment manufacturer

Quality standard 质量标准

客户的满意就是产品标准

Customer satisfaction is the product standard

PROFESSIONAL DENTAL
MANUFACTURER

MODEL
 **KS-R1**

Meet the world's first fully digital dental unit designed by the sharpest minds of the industry. Sit back and relax, your journey to the future of dentistry is just starting.

至尊待遇, 打造完美体验

Supreme treatment, creating the perfect experience

选配 - 儿童皮垫

意大利进口胎羊皮, 给儿童全方位守护

Optional Function
 Imported Italian sheepskin, protect children in all aspects.

人体工程学

我们与牙科行业的人体工程学专家合作, 设计了这款绝对出色的牙科椅, 以减轻客户的身体压力, 让他们尽可能舒适

We have collaborated with ergonomic experts in the dental industry to design this absolutely outstanding dental chair to reduce the physical pressure on customers and make them as comfortable as possible.



导轨式设计

悬挂式导轨设计配合台湾螺杆电机驱动, 超强承重能力, 椅架升降平稳, 更方便于两手操作和四手操作

The design of the suspended guide rail is matched with a Taiwan screw motor drive, which has a strong load-bearing capacity. The chair frame can be smoothly lifted and lowered, making it more convenient for both two-handed and four-handed operations.

Modern design, pleasing finishes, increased comfort, allowing patients to enjoy an unparalleled sense of comfort.

现代感的设计,令人赏心悦目的饰面,更高的舒适度,让患者享受到无比的舒适体验感



up
up
up.

选配 – 第一选位系统

减少或避免部分医疗事故的发生和纷争
为安全,做好第一道坚强防护

Optional Function
First selection system Reduce or avoid
some medical accidents and disputes,
and provide the first strong protection
for safety.



Modern design
pleasing finishes

SHOW YOUR IDENTITY AND MAKE YOU FEEL DOUBLY PRESTIGIOUS.

尽显身份, 让你倍感尊贵



先进的集成系统

前所未有的智控体验, 让治疗游刃有余
The unprecedented experience of intellectual control makes the treatment easy



意大利进口软皮

细致缝线, 高密度记忆棉, 舒适亲肤, 质感十足
Italian imported soft leather
Exquisite stitching, high-density memory foam, comfortable and skin-friendly, full of texture.



全球独创

牙科治疗椅变频直流系统
Global Innovative Dental unit Variable Frequency DC System



牙科护眼口腔灯

十孔双光可调色温感应灯
Ten bulbs double light adjustable color temperature sensor lamp



选配 - 牙科椅水气消毒系统

可选配牙科椅水气消毒系统, 适用于所有牙科椅, 水路和气路的消毒杀菌有效达到 99.999%
Optional Function
Dental chair moisture disinfection system is available. Suitable for all dental chairs. The effective disinfection of water and air tubes reached 99.999%

New simplicity and elegance aestheticism

新简雅美学主义

See simplicity as beauty, quality as luxury

MODEL
 **KS-R2**



1st BUILT-IN WIRELESS CURING LIGHT...

1S 内置无线座充光固化

独创航空铝导光棒，车规级散热方式，高照度蓝光+紫光，一秒极速固化，超轻机身，更便捷更舒适

Innovative aviation aluminum light guide rod, vehicle-level cooling system, high-intensity blue+purple light, one-second rapid curing, ultra-lightweight body, more convenient and comfortable.



有机硅超纤皮垫

耐磨防水防污，皮感柔和舒适

Organic silicon ultra-fiber cushion
 Wear-resistant, waterproof, and stain-resistant, with a soft and comfortable leather-like feel.



牙科椅水气消毒系统

可选配牙科椅水气消毒系统，适用于所有牙科椅水路和气路的消毒杀菌有效达到 99.999%

Optional Function
 Dental chair moisture disinfection system is available, Suitable for all dental chairs. The effective disinfection of water and air tubes reached 99.999%



LINKAGE CHAIR FRAME

联动补偿椅架

无与伦比的乘坐体验, 以患者舒适感出发构建牙椅科学坐感体系
 Unparalleled ride experience. The scientific sitting sensation system of dental chair was constructed based on patient comfort



X RAY DIGITAL IMAGING SYSTEM

数字化影像系统



Human-centered design philosophy creates higher value for you

人性化设计理念, 为你创造更高价值

MODEL
KS-R3



原装台湾聚醚水气管线
Original Taiwan polyether water and air tubing



灵活且宽敞

灵活的伸缩臂上, 可确保 2 手和 4 手治疗的最佳位置
只需几秒钟即可从右向左改变牙科治疗
这增加了在忙碌的一天选择治疗手法时的灵活性

The flexible telescopic arm ensures the best position for 2 - and 4-hand treatment
Change dental treatment from right to left in just a few seconds
This adds flexibility when choosing a treatment maneuver during a busy day



意大利进口电磁阀

杜绝外漏, 调节精度高, 响应快, 抗腐蚀无污染

Solenoid valve imported from Italy
Eliminate external leakage, high adjustment accuracy, fast response, corrosion resistance, energy saving and pollution-free.



选配 - 上挂式器械盘

旋转式托盘, 轻巧灵活拉杆, 让操作更便捷

Upper hanging instrument tray
Rotary tray, light and flexible pull bar, make the operation more convenient.

PRESENT THE BEAUTY OF THE ESSENCE OF RATIONALITY AND ART WITH PURE POINT, LINE, AND PLANE LANGUAGE

用纯粹的点线面语言, 呈现理性与艺术的本质之美

钢板上下平行盖
承载能力更佳, 稳固无比

Parallely covered steel plates provide improved load-bearing capacity and unparalleled stability



氛围灯工具盘

独创工具盘氛围灯, 让牙椅从此有了温度
Ambience lamp instrument tray
Innovative instrument tray ambient light, giving dental unit a warm atmosphere.



陶瓷痰盂

亚克力木纹面板, 独创直排式痰盂, 下水更加顺畅
可旋转拆卸痰盂, 方便清洁

Ceramic Spittoon
Acrylic wood-grain panel and original straight-draining spittoon make the water smoother. The spittoon can be rotated and removed for easy cleaning



旋转箱筒

可90度旋转箱筒
Side box with 90-degree rotation

人体曲线设计 缓解腰部酸痛

Ergonomic Designs
mitigating back pains

两侧环抱设计 扶托两侧肌肉

All-round Cushion Design
providing muscle support

独立脊椎支撑 分担腰椎压力

Stand-alone Spine Support
easing pressure on lower back

RATIONAL AND MINIMALIST BLACK, WHITE AND GRAY LUXURIOUS CLASSIC COMBINATION

理性质感的黑白灰, 轻奢的典范

MODEL
KS-R4 | 豪华版
 THE DELUXE EDITION



陶瓷痰盂

大理石石纹亚克力面板, 独创直排式痰盂
 下水更加顺畅。

Ceramic Spittoon
 Marble stone acrylic panel, original straight spittoon
 The water is smoother.




旋转箱筒

可 90 度旋转箱筒
 Side box with 90-degree rotation



旋转操作杆

180 度全手工具盘旋转操作杆, 更灵活更方便
 180 degrees full hand tool rotation lever



原装台湾聚醚水气管线
 Original Taiwan polyether
 water and air tubing

触摸操作面板

带九个记忆位的触摸按键控制面板, 操作按键更灵敏
 Touch operation panel. The touch
 button control panel with nine memory
 bits makes the operation buttons more
 sensitive.



丹麦力纳克电机

世界顶级品牌
 丹麦力纳克 (LINAK) 电机
 Denmark LINAK motor

可拆卸靠背

易于清洁, 内饰可以轻松拆卸以进行清洁
 为患者提供安心卫生的保障

Detachable metal backrest
 Easy to clean, the interior can be easily
 removed for cleaning. To provide patients with
 peace of mind health protection

The unique beauty of minimalistic luxury art

新包豪斯现代极简美学



MODEL
KS-R4 | 暗夜骑士
 THE DARK KNIGHT

YOUR CLINIC CAN BE UNIQUE.

Deconstructing the modern minimalist aesthetics of New Bauhaus, interpreting a unique sense of luxury artistic beauty under different lighting conditions

你的牙科诊所也可以很独特



陶瓷痰盂

亚克力面板，独创直排式痰盂，下水更加顺畅
 Acrylic panel, original straight line phlegm Meng, more smooth launching



生活需要些色彩

调色生活的画笔掌握在自己手里
 Life needs the harmony of colors, the brush of coloring life is in your own hands.



EXQUISITE LUXURY,
CAREFULLY
SELECTED PREMIUM
PRODUCTS

细致奢华, 甄选优品



MODEL
KS-R4 | 标准版
THE STANDARD EDITION

We Offer More
Than That



KS-R4 | 英伦典范
THE BRITISH STYLE

Green is the color of life, giving people a sense of vitality, inheriting the European royal design inspiration, luxury, but restrained and simple but not simple, natural warmth full of human care, the unique sense of advanced into life, its innate British royal temperament, and modern fashion elements perfect integration, creating a new era of oral medical equipment aesthetics



KS-R4 | 轻舞飞扬
THE FLY IN DANCE

According to the data, the human first time to go to the dentist nervous rate of 98%, the famous American psychology expert DR.Leonardo Ricardo and KASO Kono jointly designed the world's first caring dental integrated treatment machine, using the relationship between psychology and color, special deployment of the light dance flying series of colors, to give people peace, calm, calm. The feeling of relaxation will reduce the nervous rate of the first dental visit, especially the nervous rate of women and children, to 68.6%. KASO Kono leads the era to create a global care and humane digital oral clinic, so that your clinic will walk with love from now on.

Exquisite exterior design makes the dental chair the highlight of the clinic.

精致外观设计, 让牙椅成为诊所的亮点



MODEL
KS-R5

新型塑料外壳全新 ABS 材料, 表面 PMMA 材料层, 抗老化抗 UV, 耐化学腐蚀。人体工程学靠背, 使患者治疗过程舒适轻松

New plastic shell, new ABS material, surface PMMA material layer, anti-aging, anti-UV, anti-chemical corrosion. Ergonomic backrest for comfortable and easy patient treatment.



人体工程学靠背

靠背和座椅的组合移动支撑脊柱, 符合人体工程学的延伸, 使得患者舒适地移动到治疗位置

Ergonomic backrest
The combined movement of the backrest and the seat supports the ergonomic extension of the spine, allowing the patient to comfortably move to the treatment position.

可调节器械盘挂架

多角度调节, 可以轻松移动到正确位置, 让您拥有更佳的操作范围。

Adjustable instrument tray hanger
Multi-angle adjustment, you can easily move to the correct position. Let you have a better operating range.



可拆卸式金属靠背

可拆卸式全金属精架, 稳如磐石, 百倍安心

Removable metal backrest
The detachable all-metal chair frame is as stable as a rock and is a hundred times safer



豪华六孔感应灯

豪华 6 孔感应 LED 灯, 配上可拆卸消毒手柄

Sensor LED light
Luxury 6 bulb led lamp with sensor function. With detachable sterilization handle



可旋转陶瓷痰盂

90 度旋转方便患者吐痰, 表面光滑, 易于清洁

Rotatable ceramic phlegm
The 90-degree rotation is convenient for patients to spit, and the surface is smooth and easy to clean.



世界顶级品牌

丹麦力纳克 (LINAK) 电机
World topest brand Denmark LINAK motor



原装台湾聚醚水气管线

Original Taiwan polyether water and air tubing

FRESH AND ELEGANT STRIVING FOR EXCELLENCE

清新淡雅 精益求精

MODEL
KS-D213



A TRUSTY ORTHY CLASSIC

一个值得信赖的经典



LED dental light

Two light sources, white and yellow, can be selected to prevent the light-curing resin from curing in advance.



Removable metal backrest

The detachable all-metal chair frame is as stable as a rock and is a hundred times safer



Microfiber scratch-resistant skin

Luxurious comfort, long lasting



Removable strong and weak suction filter

Strong and weak suction flow can be adjusted



Water & Air Tubing

Original Taiwan polyether water and air tubing



Luxury curved arm

Fashionable and beautiful, the rotation is sensitive and convenient

Push-button instrument tray

pressure meter display

Professional and comfortable dental chair brings a relaxed and enjoyable experience

专业舒适的牙椅 带来轻松愉快的体验



MODEL KS-DLX303

经典系列的紧凑设计甚至适用于小型手术室, 同时也是大型牙科诊所装修的一种经济高效的方法。



控制面板采用清晰简单的按钮布局、图形和高品质的表面处理, 我们保持简单, 因此您可以完全专注于您的工作, 不会出现任何干扰或错误

The multifunction expansion table is located on a flexible telescopic arm to ensure optimal placement for 2- and 4-hand treatments. Change dental treatment from right to left in just a few seconds.



让你的病人得到便捷

90度旋转方便患者漱口, 表面光滑, 易于清洁

The 90 degree rotation is convenient for the patient to gargle. Smooth surface, easy to clean.



多功能脚踏

方便控制椅位靠背升降、俯仰, 一键冲痰等功能

Multifunction foot pedal
Convenient control of the chair backrest lifting, pitching, one-key phlegm flushing and other functions.



丹麦力纳克电机

世界顶级品牌: 丹麦力纳克 (LINAK) 电机

The world's top brand
Denmark LINAK motor



原装台湾聚醚水管线

Original Taiwan polyether water and air tubing

Personalized adjustment function meets the needs of different groups of people

个性化调节功能 满足不同人群的需求

MODEL
KS-DLX301



A TRUSTY ORTHY CLASSIC

一个值得信赖的经典



豪华弯型转臂

时尚美观, 旋转灵敏便捷

Luxury curved arm
Fashionable and beautiful, the rotation is sensitive and convenient

触摸按键器械盘 (带气压表)

Touch button instrument tray (pressure meter display)



靠背控制系统

便于医生和病人操作的靠背控制系统
专利号: 201721249324.4

Backrest control system
Backrest button system convenient for doctors and patients to operate Patent No. 201721249324.4



机箱线路优化

水电分离, 排列整齐, 易于维修


Side box circuit optimization
Separation of water and electricity, neatly arranged easy to maintain.



辅控按键系统

方便助手对牙椅常用功能控制

Assistant Control system
Convenient for assistants to use dental chair functions.



原装台湾聚醚水气管线

Original Taiwan polyether water and air tubing

SIMPLE DESIGNS CAN ALSO MEET THE REQUIREMENTS.

简易的设计亦能满足需求

MODEL
KS-D108



选配 - 靠背控制系统

便于医生和病人操作的靠背控制系统
专利号 :201721249324.4


Optional
Backrest control system convenient for dentist
Patent No.201721249324.4



简洁机箱

小巧的机箱为医师、护士提供良好的观察视角及操作空间

Simple side box
The compact chassis provides doctors and nurses with a good viewing angle and operating space.



原装台湾聚醚水管线

Original Taiwan polyether water and air tubing

高效的工作空间

在设计时充分考虑了您和您的团队的需求
由于设计中融入了人体工学原理, 工作变得非常简单和直观工作日对身体的负担也被最小化

Fully considering the needs of you and your team during design
Due to the integration of ergonomic principles into the design, work has become very simple and intuitive, and the burden on the body during workdays has been minimized

UNIQUE STREAMLINED APPEARANCE

独特流线型外观

Incorporating modern streamlined design elements, the overall shape of the dental chair is not only practical but also adds a sense of fashion and taste to the entire clinic, making it an integral part of the environment. This enhances the patients' visiting experience.



MODEL
KS-R6

全新设计 灵活机动 便捷安装 省时省力



创新活动座椅和机箱旋转结构，无需拆卸，让安装和移动变得前所未有的简单，轻松解决传统牙椅进窄门的困扰，不论是狭窄走廊还是小型诊室，轻松适应各种空间，减少拆卸安装所需的时间和劳力，为牙科诊所提供无限方便。

让我们的牙椅成为您牙科诊疗的得力助手！

The innovative activity seat and cabinet rotation structure eliminate the need for disassembly, making installation and relocation unprecedentedly simple. This solution effectively addresses the challenges of maneuvering traditional dental chairs through narrow doors. Whether in tight corridors or compact clinics, our chair easily adapts to various spaces, significantly reducing the time and labor required for setup.

Let our dental chair be your reliable assistant in dental practice!



INNOVATIVE DESIGN WITH FLEXIBILITY
CONVENIENT INSTALLATION - SAVES TIME AND EFFORT

MODEL
KS-R6



" HUMANIZED
COMFORT
EXPERIENCE

"FLEXIBLE
SPACE"



轻松折叠, 省力省时, 创造更多灵活空间

"Easy to fold,
saving effort
and time,
creating
more flexible
space."

OPTIONS FOR DENTAL CHAIR

牙科治疗机可选配件



KS-TM111 可升降小推车
Adjustable Treatment Trolley



牙科椅水气消毒系统
Dental Chair Moisture
Disinfection System



KS-W111B 美白仪
Teeth Whitening Accelerator



内置光固化
Built-In Light Cure



内置洁牙机
Built-In Scaler



KS-CM109 数字化口腔观察仪
Digital Dental Camera



高低速手机套装
Handpiece Set



KS-AirJet 1 喷砂洁牙手机
Powder Jet Handpiece



内置电动马达
Built-In Electric Motor



KS-Q8 1秒光固化
One Second Light Cure



KS-ED07 根管治疗仪
Wireless Endo Motor



KS-IS4100 种植机
Dental Implant Unit



KS-OT03 操作手推车
Operating Trolley



KS-W117A 立式美白仪
Teeth Whitening Machine



KS-W119 立式美白仪
Teeth Whitening Machine



KS-R7 根管测量仪
Apex Locator



KS-A1008 根管测量仪
Apex Locator



KS-E1002 二合一根管治疗仪
Endo Motor With Apex Locator



KS-S7 超声波洁牙机
Ultrasonic Scaler



KS-EM03 电动马达
Dental Electrical Motor



KS-TYPE 2 超声波洁牙机
Ultrasonic Scaler

OILFREE AIR COMPRESSOR

无油静音空压机

- 无油科技舒适静音: 选用优质铁壳消音器, 静音效果更出色
- 无油维护: 无需加油润滑, 气源洁净, 有效提升用户设备的使用效果和使用寿
- 全铜加长电机: 动力强劲、经久耐用、高性能、低能耗
- 大口径散热孔: 发热缓慢、风扇快速散热、可长时间开机使用
- 配备安全: 内置过热保护器, 加强电机保护
- 进口皮碗: 更耐磨, 更静音
- 售后保障: 整机质保1年, 免除用户后顾之忧



| | KS-C1001 | KS-C1002 | KS-C1003 | KS-C1004 | KS-C1005 | KS-C1006 |
|-------------------------------|-----------------------|-------------------------|------------------------|-------------------------|-------------------------|-------------------------|
| 电 源 Voltage(V) | 110v(60Hz)/220v(50Hz) | 110v (60Hz) /220v(50Hz) | 110v (60Hz)/220v(50Hz) | 110v (60Hz) /220v(50Hz) | 110v (60Hz) /220v(50Hz) | 110v (60Hz) /220v(50Hz) |
| 电 流 Current(A) | 2.8 | 7.6/3.8 | 5.6 | 15.2/7.6 | 8.4 | 11.5 |
| 功 率 Power(Kw) | 600 | 850 | 1200 | 1700 | 1800 | 2550 |
| 公称容积流量 Air flow(L/min) | 90 | 110 | 180 | 220 | 270 | 330 |
| 额定排气压力 Capacity Pressure(Mpa) | 0.8 | 0.8 | 0.8 | 0.8 | 0.8 | 0.8 |
| 储气罐容积 Tank (L) | 35 | 35 | 60 | 60 | 90 | 90 |
| 噪 音 Noise(dB) | 63 | 65 | 68 | 70 | 72 | 72 |
| 净 重 Net weight(KG) | 22 | 25 | 40 | 45 | 60 | 65 |
| 毛 重 Gross weight (KG) | 24 | 27 | 55 | 60 | 80 | 80 |
| 机 身 尺 寸 Unit Dimension(mm) | 410x410x600 | 410x410x620 | 800x370x690 | 800x370x710 | 950x370x690 | 950x370x710 |
| 包 装 尺 寸 Packing Dimension(mm) | 450x450x660 | 450x450x660 | 990x520x820 | 990x520x820 | 1140x520x820 | 1140x520x820 |
| 连接牙椅数量 For Dental Units | 1 | 1~2 | 2~3 | 3~4 | 4~5 | 5~6 |

SUCTION PUMP

牙科抽吸机



| | SPL02 | SPL05 |
|-------------------------|------------------|------------------|
| Model | SPL02 | SPL05 |
| Supported Dental Chairs | 2~3 | 5~8 |
| Input Volt | 220V(1~)/50~60HZ | 220V(1~)/50~60HZ |
| Power | 750W | 1850W |
| Max Current | 5.5A | 7.5A |
| Vacuum | -15KPA | -21KPA |
| Maximum Airflow | 1100L/MIN | 2000L/MIN |
| Frequency Conversion | NO | YES |
| Noise Level | 53DB(A) | 55DB(A) |
| Weights | 23KG | 45KG |
| Dimension (LxBxH) | 286x500x600 MM | 488x400x790 MM |

DENTAL CENTRALIZED SUCTION SYSTEM

中央负压系统



| | SPL12 |
|-------------------------|------------------|
| Model | SPL12 |
| No. of vacuum pumps | 4 |
| Supported Dental Chairs | 40~50 |
| Input Volt | 380V/50HZ (3~) |
| Power | 15.5 |
| Control method | PLC CONTROL |
| Frequency Conversion | YES |
| Vacuum | -24 |
| Maximum Airflow | 20000 |
| Noise Level | 70 |
| Weights | 1325x650x1370 MM |
| Dimension (LxBxH) | 850x700x1570 MM |

KS-KU AUTOCLAVE

高压灭菌器



| Technical specifications | 18L | 23L |
|-----------------------------|--|-------------|
| Chamber (mm) | φ247 x 350 | φ247 x 450 |
| Overall dimensions (W*H*D) | 490*455*600 | 490*455*690 |
| Net Weight (kg) | 47 | 53 |
| Nominal power (VA) | 1750 | 1750 |
| Rated Voltage | 230~240V 50 Hz | |
| Sterilization temperatures | 121° C /134° C | |
| Capacity of the distilled | 2.5 L (Water at level Max.) Approx. 0.5 L (Water at level Min.) Approx. | |
| Circuit breaker | F16A /400 V | |
| Operation temperature | 5° C ~ 40° C | |
| Operation relative humidity | Max. 80%, non-condensing | |
| Max. Noise level | <70 dB | |
| Atmospheric pressure | 76 kPa ~ 106 kPa | |

SEALING MACHINE

KS-F-III 封口机



| | |
|----------------------|---------------------|
| Voltage | AC 110~230V 50~60Hz |
| Dimensions | 370x440x 195mm |
| Packaging dimensions | 430x340x250mm |
| Sealing width | 0~250mm (max) |
| Sealing temperature | 220° C±10° C |
| Total power | 500W |
| Indentation width | 10mm |
| Net weight | 5.5Kg |
| Fuse tube | T3A |

AUTOCLAVE

高压灭菌器



| Model | KS-BE18L | KS-BE23L |
|---------------------------|---------------------|---------------|
| Voltage | 110V/240V,50Hz/60Hz | |
| Power | 1500W | 2000W |
| Volume | 18L | 24L |
| Chamber Dimension | Ø245x355mm | Ø245x465mm |
| Relative Humidity range | << 85% | |
| Ambient temperature range | +5° C ~ +40° C | |
| Machine Dimension | 590x415x410mm | 590x415x410mm |
| Weight | 42Kg | 45Kg |

KS-WATER DISTILLER

KS-BV-1 蒸馏水机



| | |
|--|--------------------|
| Rated voltage | 220~240VAC 50Hz |
| Rated power | 750W |
| Inner cavity size | Diameter 180x200mm |
| Overall dimensions (length x width x height) | 290x290x390mm |
| Packaging dimensions (length x width x height) | 300x300x500mm |

X-RAY MACHINE

High resolution, providing clearer images, adaptable to different examination needs, and generally more stable.

KS-X105 / KS-X106 X-ray Machine

| | | |
|------|---------------------|------------------|
| 辐射焦距 | Tube Focus | 1.5 |
| 球头型号 | Tube Model | XD2-1.4/85 |
| 球头电压 | Tube Voltage | 60KV |
| 球头电流 | Tube Current | 8mA |
| 固定滤波 | Inherent Filtration | 0.75mmAL |
| 辐射泄漏 | Emissivity Leakage | 1 metres0.1mGy/h |
| 曝光时间 | Exposure Time | 0.2-3.2S |



AT-301

高清口腔数字传感器
Digital intraoral imaging system

iSensor(iTech Sensor)

高清口腔数字传感器
Digital intraoral imaging system Sensor


KUNLUN

一台拍的更高清的
智能CBCT影像设备



TOSHIBA HIGH-END BALL TUBE

The use of CBCT applications more high-end ball tube model, continuous shooting performance is always stable, the best partner for high-definition images



SUPPORTS FOV UP TO 20X12 FIELD OF VIEW

It can provide 14x10, 15x9, 15x10.5, 16x1120*12 multi-size FOV to meet the clinical needs of small, medium and large field of vision in different complex surgical scenarios



EXTRA-LARGE FLAT PANEL DETECTOR

Amorphous silicon ultra-large size, non-splicing imaging, shooting image stability, anti-attenuation, long life, high-definition image quality



INDEPENDENT RESEARCH SOURCE SYSTEM

Intelligent real-time temperature control, overtemperature warning protection and air bag design, to avoid the deformation of the ball tank, durable as new

DENTAL INSTANT CAMERA

The body is small, light weight and nearly no radiation. It has excellent image quality, portable storage, save more space.



Vision X-A



Vision X-B

X-ray Unit Specification

| | |
|----------------|----------------------|
| Tube Voltage | 60(70) kV |
| Tube Current | 2.0 (1.0,1.5,2.5) mA |
| Exposure Time | 0.1-2s |
| Tube Focus | 0.8x0.8(0.4x0.4)mm |
| Battery | DC 14.8V,1300mA |
| Rated Power | 120(100-175) VA |
| Charger Input | AC100v-240V |
| Charger Output | DC 16.8 V |

Digital Sensor Specification

| Model | COX-01 | COX-02 |
|--------------------------|-----------|-----------|
| Effective pixel (pixel) | 1200x1600 | 1402x1888 |
| Effective area size (mm) | 22.2x29.6 | 26x35 |
| Resolution ratio (lp/mm) | >10 | >10 |
| Output | USB 2.0 | USB 2.0 |
| Power (mW) | 350 | 350 |
| ESD (kV) | 8 | 8 |
| Service temperature (°C) | 0~40 | 0~40 |

PORTABLE X-RAY UNIT

Due to its small size, it is portable and easy to carry. It can also operate in narrow spaces, making it very flexible. Additionally, it has a fast imaging speed.



KS-PX01



KS-PX02



KS-PX03



KS-PX05



KS-PX06



KS-PX07

| 球管焦点 | Tube Focus | 0.8 | 0.4 | 0.4 | 0.8 | 0.8 | 0.8 |
|------|---------------------|---------------|---------------|---------------|-----------|-----------|-----------|
| 球管电压 | Tube Voltage | 60KV | 60KV | 60KV | 60KV | 60KV | 60KV |
| 球管电流 | Tube Current | 1.5mA | 1.5mA | 2mA | 2mA | 2mA | 2mA |
| 电池容量 | Battery Capacity | 5000mA | 3500mA | 900mA | 1600mA | 1600mA | 1600mA |
| 固有过滤 | Inherent Filtration | 0.5mmAL | 1mmAL | ≥ 0.75mmAL | 0.5mmAL | 0.5mmAL | 0.5mmAL |
| 辐射泄漏 | Emissivity Leakage | 1m ≤ 0.1mGy/h | 1m ≤ 0.1mGy/h | 1m ≤ 0.1mGy/h | 0.05mGy/h | 0.05mGy/h | 0.05mGy/h |
| 曝光时间 | Time Of Exposure | 0.01-2.5S | 0.01-2.5s | 0.1-2.0S | 0.1-2S | 0.1-2S | 0.1-2S |
| 频率 | Frequency | 23KHz | 23KHz | 80KHz | 30kHz | 30kHz | 30kHz |

BUILT-IN WIRELESS CURING LIGHT

Experience Ultimate Convenience with the All-New Built-in Wireless Technology

1S

ONE SECOND FASTCURING

SmartCure+



Aluminum alloy ultra-small elbow

- 1) Aviation aluminum alloy, strong durability;
- 2) Diameter is only 8mm, 45-degree bend angle, easy to insert into the mouth;
- 3) All materials are high-temperature resistant;
- 4) Five-pin structure, better and more stable contact.



wireless

1 SECOND FAST CURING

Aluminum alloy detachable ultra-small elbow

360-degree rotation

Time display

5 working modes



Technica Specifications

| | | |
|------|-----------------------|-----------------------------|
| 输入电压 | Charger Input Voltage | DC5V |
| 电池电压 | Battery Voltage | 3.7-4.2V |
| 电池容量 | Battery Capacity | 300mAh |
| 波长 | Wavelength | 400-460nm |
| 照明范围 | Illuminance | 1000-4000mW/cm ² |

Five working modes


- Normal Mode
- Fixed second mode
- Pulse Mode
- Gradual Mode
- Caries Mode

3D Intraoral Scanner

Digitize your dentistry with us

ScanMax

Bid adieu to traditional dental impressions. Experience the future of dental diagnostics



A.I.

A.I. Scan: A.I. scanning will automatically filter out unnecessary data, like tongue and lips during scanning, which can keep the scanning process smooth.

Free Scan

Real-time global alignment
No scan path dependence

Scan Depth 25MM

Higher scan efficiency
Stronger inlay and implant scanning capabilities

Safety

Gesture control - reduce the risk of cross infection while improving Scan workflow efficiency

高精度
高精度是数字牙科工作流程的基础
Accuracy: High accuracy is the basis of the digital dentistry workflow.

< 20_{μm}

捷 - 极速体验
全牙弓扫描不到 40 秒
Speed: Provides you the smooth scanning experience.

< 30_s

触屏电脑一体机配置 (选配 optional)
Touchscreen all-in-one computer



| | |
|----------------|-------------------------|
| 操作系统 OS | Windows10 (64 位) |
| 处理器 CPU | 第 11 代 intel i5 及以上 |
| 内存 Memory | 16GB |
| 硬盘 Hard Disc | 512G |
| 显卡 GPU | NVIDIA GTX 3060 6GB 及以上 |
| 显示器 Resolution | 1920x1080, 60Hz |
| 接口 USB Port | 至少 2 个 USB3.0 |

技术参数
Specification

| | |
|----------------------|----------------------------|
| 产品名称 Product Name | 口内扫描仪 3D Intraoral Scanner |
| 文件导出类型 Output Format | STL, PLY, OBJ |
| 精度 Accuracy | 20um |
| 景深 Scan Depth | 25mm |
| 电源 Power supply | DC12V 2A |
| 电缆长度 precision | 1.8m |
| 重量 Weight | 270g |
| 尺寸 Size | 264x37x46mm |
| 扫描头 Tip Size | 20x17mm |
| 光源 Light Source | LED |
| 手势控制 Gesture control | 是 YES |
| AI 技术 A.I. Tech | 是 YES |

数字化口内智能扫描牙椅
KASO digital 3D scanning chair



ORAL TREATMENT MICROSCOPE MICRO-X



Built-in long-travel zoom objective lens
F=250mm

A large zoom objective lens suitable for surgery, dentistry, and orthopedics. The focusing range covers the entire depth of the oral cavity. The large working distance achieved by the large zoom makes it easy to adjust the zoom focus.

Excellent optical system

The self-developed imaging optical system is adopted, demonstrating the perfect combination of optical design and digital imaging technology. It has high definition, strong three-dimensional sense, and large depth of field. It can achieve 180° naked-eye operation, and the intraoperative hand-eye synchronization rate reaches 99.8%, 5 times compared with the magnification of the traditional microscope, it breaks the traditional method, frees the doctor's eyes, satisfies the observation of different tooth positions, and improves the doctor's work efficiency.



-  **White light mode**
Using imported medical LED light source, high color rendering index, the color is closest to reality, and the illumination spot has good uniformity.
-  **Orange light mode**
Used for filling light-curing resin materials to prevent the resin materials from curing too quickly.
-  **Green light mode**
See the tiny nerves and blood vessels clearly in the surgical blood environment, and improve the contrast of blood vessels and bleeding points under blood flow conditions.



1 Pendulum system
The mirror body tilts 45° left and right. When you observe from different positions, the field of view remains unchanged.



2 Handle
Ergonomically designed handle, easy to operate.



3 Touch display
The fuselage has a built-in touch screen design, and an external mouse can also be connected to set the system's menu function parameters and other focus and focus settings.




4 left control handle
The main operating functions include flipping the microscope up and down, left and right at four different angles, viewing preview and return keys to meet the needs of treatment from different angles.



5 Right control handle
The main operating functions are manual focus, one-button focus and auto-focus, ZOOM continuous zoom system adjustment, defogging and lighting and other key functions.



6 Large zoom objective
A large zoom objective lens suitable for surgery, dentistry, and orthopedics. The focusing range covers the entire depth of the oral cavity. The large working distance achieved by the large zoom makes it easy to adjust the zoom focus.



7 The mirror body is designed with a large tilt angle in the front and back.
The mirror body can be swung forward at 45° and back at 90°, which can adapt to various special treatment positions and working angles that are difficult to observe.



8 4K imaging module
Collect 4K ultra-high-definition photos and images and output them to a 4K monitor. The exquisite images will make you more comfortable in academic exchanges and doctor-patient communication.



180° wide field of view and large depth of field

There are no position requirements or eyepieces to constrain the eyes, and communication between surgeons and nurses on the same stage is more accurate. Breaking the traditional method, freeing up the doctor's eyes and operating space, allowing for more flexible surgery. Ultra-large magnification and ultra-large depth of field bring more refined treatment and reduce surgical trauma.

XM110 standard version

| | |
|-------------------------------|---------------------------------|
| Magnification factor | 1x~3x |
| Overall | 3.0x~14x |
| Magnification and zoom system | Continuous zoom, up to 77 times |

"MINGJING" SERIES DENTAL LAMPS

Caring for the health of dental doctors starts with a lamp


“明镜”系列口腔灯——关爱口腔医生健康从一盏灯开始



Ten bulbs double light adjustable color temperature sensor lamp



Six bulbs double light adjustable color temperature sensor lamp



All use original imported LED light sources with "EyeComfort" technology. No visible flicker, reliable quality and durable.



RGO Exempt Blu-ray International standard, strength eye protection



High color temperature consistency Uniform luminous effect



LCD touch control design intuitive and clear Simple and easy to use, it provides doctors and nurses with a convenient operating experience.

| | | |
|--------------------|---------------------------|--------------------------|
| Model | 6 bulbs | 10 bulbs |
| Input voltage | AC/DC 12V/24V | |
| Power | 7W | 10W |
| Frequency | 50/60 HZ | |
| Illumination | 10000-30000 (adjustable) | 10000-50000 (adjustable) |
| Colour temperature | 3000-5700K(adjustable) | |
| CRI | Ra95 | |
| Light spot size | 7*140 MM (distance 70 MM) | |



全光谱



93** CRI(Ra) 高显色指数

Adopting unique "FREESYNC" light source synchronized true color technology According to laboratory testing, the color rendering index of KASO LED dental lights exceeds 93. It improves tooth color and provides dental experts with clearer detailed judgments.



DOCTOR CHAIR

医生椅

Everyone is different Every patient, clinic, team So we put our products in front of you, not the other way around. There's a way. There's a way. Your way



KS-004



KS-007



KS-008



KS-016



KS-026 (Plastic)



KS-026



KS-029



KS-034



KS-041



KS-042



KS-045



KS-055

MICROFIBER SCRATCH-RESISTANT SKIN

超纤耐刮皮

Scratch-resistant leather color (Model:R3,R4)



070-04



070-12



070-13



070-16



070-17



070-18



070-19



070-20



070-21



070-22



070-23



070-24



070-25



070-26



070-29



070-31



070-32



070-33



070-36



070-38



070-41



070-45

OPTIONAL COLOR FOR

PU 皮套

PU Cushion color (Model:213.106.108.301.303.R5)



P1



P2



P3



P4



P5



P6



P7



P8



P9



P10



P11



P12



P13



P14



P15



P16



P17



P18



P19



P20



P21



P22



P23