

FG 高速车针 High Speed Carbide Burs



FG 球形 ROUND HEAD(SPHERICAL)											
US NO.	ISO NO.	1/4	1/2	1	2	3	4	5	6	7	8
FG	500314001001	005	006	008	010	012	014	016	018	021	023



FG 倒截锥形 INVERTED CONE HEAD (INVERTED.TRUNCATED CONICAL)										
US NO.	ISO NO.	33 1/2	34	35	36	37	38	39	40	41
FG	500314010001	006	008	010	012	014	016	018	021	023



FG 梨形 & 加长型 PEAR HEAD (HEMISPHERICAL.INVERTED CONICAL) & LONG HEAD										
US NO.	ISO NO.	329	330	330L	331	331L	332	332L	333	333L
FG	500314237001	006	008	—	010	—	012	—	014	—
FGL	500315238001	—	—	008	—	010	—	012	—	014



FG 圆柱裂沟形 & 加长型 STRAIGHT FISSURE HEAD (CYLINDRICAL) & LONG HEAD											
US NO.	ISO NO.	54	55	56	56L	57	57L	58	58L	59	59L
FG	500314107006	005	006	008	—	010	—	012	—	014	—
FGL	500316109006	—	—	—	008	—	010	—	012	—	016



FG 圆柱裂沟形(带横切) & 加长型 STRAIGHT FISSURE HEAD (CYLINDRICAL) & LONG HEAD											
US NO.	ISO NO.	556	556L	557	557L	558	558L	559	559L	560	560L
FG	500314107007	008	—	010	—	012	—	014	—	016	—
FGL	500315109007	—	008	—	010	—	012	—	014	—	016



FG 圆柱半球裂沟形 STRAIGHT FISSURE HEAD WITH ROUNDED END (HEMISPHERICAL.CYLINDRICAL)					
US NO.	ISO NO.	1156	1157	1158	1159
FG	500314137006	008	010	012	014



FG 圆柱半球裂沟形 STRAIGHT FISSURE HEAD WITH ROUNDED END (HEMISPHERICAL.CYLINDRICAL)			
US NO.	ISO NO.	244	245
FG	500314233006	006	010

FG 高速车针 High Speed Carbide Burs



FG 圆柱半球裂沟形(带横切) & 加长型 STRAIGHT FISSURE HEAD WITH ROUNDED END (HEMISPHERICAL.CYLINDRICAL) & LONG HEAD										
US NO.	ISO NO.	1556	1557	1557L	1558	1558L	1559	1559L	1560	1560L
FG	500314137007	008	010	—	012	—	014	—	016	—
FGL	500315139007	—	—	010	—	012	—	014	—	016



FG 截锥裂沟形 & 加长型 TAPERED FISSURE HEAD & LONG HEAD										
US NO.	ISO NO.	168	169	169L	170	170L	171	171L	172	172L
FG	500314168006	008	009	—	010	—	012	—	016	—
FGL	500315170006	—	—	009	—	010	—	012	—	016



FG 截锥裂沟形 & 加长型(带横切) TAPERED FISSURE HEAD (TRUNCATED CONICAL) & LONG HEAD										
US NO.	ISO NO.	699	699L	700	700L	701	701L	702	702L	703
FG	500314168007	009	—	010	—	012	—	016	—	021
FGL	500315170007	—	009	—	010	—	012	—	016	—

临床

Crown Prep.	FG 1957	FG 1958
Fine . Cross	137/010	137/012

Scaling Bur 去结石车针	Perio
	009

Scaling Bur 去结石车针	Flame
	012

Endo End Bur 安全端车针	Endo-Z	Endo-Z	Endo-Z	Endo-Z
	016	016	016	016

Extraction Bur 拔牙车针	Zekrya	Zekrya	Zekrya	Zekrya

RA 低速车针

Low Speed Carbide Burs



RA 球形 ROUND HEAD(SPHERICAL)

US NO.	ISO NO.	1/2	1	2	3	4	5	6	7	8
RA	500204001001	006	008	010	012	014	016	018	021	023



RA 倒截锥形 INVERTED CONE HEAD (INVERTED,TRUNCATED CONICAL)

US NO.	ISO NO.	33 1/2	34	35	36	37	38	39	40	41
RA	500204010001	006	008	010	012	014	016	018	021	023



RA 梨形 & 加长型 PEAR HEAD (HEMISPHERICAL,INVERTED CONICAL) & LONG HEAD

US NO.	ISO NO.	329	330	330L	331	331L	332	332L	333	333L
RA	500204237001	—	008	—	010	—	012	—	—	—
RAL	500315238001	—	—	—	010	—	012	—	—	014



RA 圆柱裂沟形 & 加长型 STRAIGHT FISSURE HEAD (CYLINDRICAL) & LONG HEAD

US NO.	ISO NO.	55	56	56L	57	57L	58	58L	59	59L	60	60L	61	62	63
RA	500204107006	006	008	—	010	—	012	—	014	—	016	—	018	021	023
RAL	500204109006	—	—	009	—	010	—	012	—	014	—	016	—	—	—



RA 圆柱裂沟形(带横切) & 加长型 STRAIGHT FISSURE HEAD (CYLINDRICAL) & LONG HEAD

US NO.	ISO NO.	556	556L	557	557L	558	558L	559	559L	560	560L	561	562	563
RA	500204107006	008	—	010	—	012	—	014	—	016	—	018	021	023
RAL	500315109007	—	008	—	010	—	012	—	014	—	016	—	—	—

RA 低速车针

Low Speed Carbide Burs



RA 截锥裂沟形(带横切) & 加长型 TAPERED FISSURE HEAD (TRUNCATED CONICAL) & LONG HEAD

US NO.	ISO NO.	699	700	700L	701	701L	702	703	704
RA	500204168007	009	010	—	012	—	016	021	023
RAL	500204170007	—	—	010	—	012	—	—	—



RA 圆柱半球裂沟形(带横切) & 加长型 STRAIGHT FISSURE HEAD WITH ROUNDED END (HEMISPHERICAL,CYLINDRICAL) & LONG HEAD

US NO.	ISO NO.	1556	1557	1557L	1558	1558L	1559	1559L	1560	1560L
RA	500204137007	008	010	—	012	—	014	—	016	—
RAL	500204139007	—	—	010	—	012	—	014	—	016

RA 碳钢低速车针

Stainless Steel Low Speed Burs



RA 球形 ROUND HEAD(SPHERICAL)


US NO.	ISO NO.	2	3	4	5	6	7	8
RA		010	012	014	016	018	021	023


种植用球钻 快速处理软组织 RA XL 30.0mm Implant use

Ball	RA XL D10	RA XL D20	RA XL D30	RA XL D40	RA XL D50
	1.0	2.0	3.0	4.0	5.0

骨科用球钻 HP XL 70.0mm Bone Surgical

Ball	HP XL D10	HP XL D20	HP XL D40
	1.0	2.0	4.0



车针柄型号 Specification	车针柄直径 Diameter	车针柄长度 Length
FG(Friction Grip) 	1.60mm	19.0mm
SS(FG Short Shank) 	1.60mm	16.0mm
FG SL(FG Surgical Length)  (外科长度)	1.60mm	25.0mm
RA(Right Angle) 	2.35mm	22.5mm
RA SL(RA Surgical Length)  (外科长度)	2.35mm	25.0mm
HP(Handpiece) 	2.35mm	44.5mm

HP 技工用车针 Handpiece Carbide Burs (Laboratory Use)



HP 球形 ROUND HEAD(SPHERICAL)													
US NO.	ISO NO.	1/2	1	2	3	4	5	6	7	8	9	10	11
HP	500104001001	006	008	010	012	014	016	018	021	023	025	027	31



HP 倒截锥形 INVERTED CONE HEAD (INVERTED, TRUNCATED CONICAL)									
US NO.	ISO NO.	34	35	36	37	38	39	40	41
HP	500104010001	008	010	012	014	016	018	021	023



HP 圆柱裂沟形 STRAIGHT FISSURE HEAD (CYLINDRICAL)									
US NO.	ISO NO.	56	57	58	59	60	61	62	63
HP	500104107006	008	010	012	014	016	018	021	023



HP 圆柱裂沟形(带横切) & 加长型 STRAIGHT FISSURE HEAD (CYLINDRICAL) & LONG HEAD															
US NO.	ISO NO.	556	556L	557	557L	558	558L	559	559L	560	560L	561	562	563	
HP	500104107007	008	-	010	-	012	-	014	-	016	-	018	021	023	
HPL	500104110007	-	008	-	010	-	012	-	014	-	016	-	-	-	



HP 截锥裂沟形(带横切) & 加长型 TAPERED FISSURE HEAD (TRUNCATED CONICAL) & LONG HEAD												
US NO.	ISO NO.	699	699L	700	700L	701	701L	702	702L	703	703L	704
HP	500104168007	009	-	010	-	012	-	016	-	021	-	023
HPL	500104171007	-	009	-	010	-	012	-	016	-	021	-